



2022

Eco-Suture Catalog

aspen
veterinary
resources





Committed to Quality.

Dedicated to Animal Health.

For almost 30 years, Aspen has provided our valued customers with superior products formulated to optimize animal health.

At Aspen, our products are designed specifically for animal health professionals to ensure the animals they treat receive the very best medical attention possible. Every product – from wound closure, to disease prevention, to reproduction and nutrition – has been thoroughly researched and hand-picked for our customers and the important clients they serve.

We are committed to expanding our portfolio using only the highest-quality materials and cost-effective solutions. That's our promise.

Thank you for trusting Aspen with the health of the animals in your care.



Table of Contents

THE ASPEN DIFFERENCE – UNPARALLELED WOUND CARE	4
ASPEN ECO-SUTURE PACKAGING	6
ECO-SUTURE OVERVIEW	8
ECO-CRYL™ PGLA	9
ECO-MAX™ PDO	10
ECO-LON™ NYLON	11
ECO-WEB™ PGCL	12
COMMON USES	14



The Aspen Difference - Unparalleled Wound Care

At Aspen, we understand that choosing the right suture for each procedure is as unique as the veterinarians using them. Designed specifically for veterinarians, Aspen offers a diverse selection of swaged sutures. Our absorbable and non-absorbable suture profiles provide wound closure options that reduce the risk of infection and promote healing.

Exceptional Quality

Aspen Eco-Sutures come with proprietary 300 series stainless steel needles, making them more durable and less brittle than other suture products. Our needles are designed to meet the demands of procedures involving different tissue types.


Purposeful Packaging

Aspen Eco-Sutures are sold in economical boxes of 12 to reduce waste and manage inventory costs.

Made in the USA

Aspen Eco-Sutures are proudly manufactured in the USA, in an environmentally conscious FDA-registered facility with class 10,000 ISO-7 clean rooms to ensure sterility and product shelf life.






ASPEN ECO-SUTURE PACKAGING

PACKAGING OVERVIEW

Every aspect of Aspen's Eco-Suture line was designed with the end user – our valued customers – in mind. Our suture boxes are aligned with industry color coding for easy identification and each box clearly identifies the contents inside.


These sutures stand up to Aspen's brand values – to provide cost-effective solutions without compromising quality.



COLOR CODED PACKAGING
Matching the industry standards.

ASPEN ECO-SUTURE PACKAGING

PACKAGING FEATURES



Suture Box

- Industry color coded for easy identification
- Easy to open box with over-wrap pull tab for effortless opening
- Made from recyclable materials

Suture Peel Pouch


- Foil pouch offers complete moisture and UV barrier
- Easy to open peel design
- Individual foil pouch to reduce waste and steps to get product onto a sterile field

Protective Foil Pouch



EASY TO USE

Aspen Eco-Sutures are as easy to use as 1, 2, 3.



ASPEN ECO-SUTURE OVERVIEW

Aspen Eco-Suture Product Name

ECO-CRYL™ PGLA

- Coated
- Braided
- Polyglactin 910 (PGLA)
- Synthetic Absorbable Suture
- Suture color: Violet

All Species



Material Comparison

Covidien | Polysorb[^]
Ethicon | Vicryl*
Covetrus | PGA^{^^}
VetOne | One-Sorb^{**}
VPL | Polymend ES^{^^^}

ECO-MAX™ PDO

- Monofilament
- Polydioxanone (PDO)
- Synthetic Absorbable Suture
- Suture color: Violet

All Species




Covidien | Maxon[^]
Ethicon | PDS II*
Covetrus | Monosorb^{^^}
VetOne | One-Dox^{**}
VPL | Monomend Max^{^^^}

ECO-LON™ NYLON

- Monofilament Nylon
- Non-Absorbable Suture
- Suture color: Blue

All Species



Covidien | Monosof[^]
Ethicon | Ethilon*
Covetrus | Monofil^{^^}
VetOne | One-Lon^{**}
VPL | NY-STA^{^^^}
















ECO-WEB™ PGCL

- Monofilament
- Poliglecaprone 25 (PGCL)
- Synthetic Absorbable Suture
- Suture color: Violet








All Species









Covidien | Caprosyn/Biosyn[^]
Ethicon | Monocryl*
Covetrus | Sorbocryl^{^^}
VetOne | One-Cryl^{**}
VPL | Monomend MT^{^^^}

Description	Item Code	Needle Type	Material	Suture Color	Size	Length	Mfg Code	Quantity
ECO-CRYL™ PGLA								
CP-1 36 mm ½ circle 	21293334	 Reverse Cutting	Braided	Violet	0	30 in (75 cm)	AB-V467	12s
CT-1 36 mm ½ circle 	21293324	 Taper Point	Braided	Violet	2-0	30 in (75 cm)	AB-V339	12s
CT-2 26 mm ½ circle 	21293322	 Taper Point	Braided	Violet	2-0	30 in (75 cm)	AB-V333	12s
	21293321	 Taper Point	Braided	Violet	3-0	30 in (75 cm)	AB-V332	12s
FSL 30 mm ¾ circle 	21293338	 Reverse Cutting	Braided	Violet	0	30 in (75 cm)	AB-V587	12s
FS-1 24 mm ¾ circle 	21293330	 Reverse Cutting	Braided	Violet	2-0	30 in (75 cm)	AB-V453	12s
	21293329	 Reverse Cutting	Braided	Violet	3-0	30 in (75 cm)	AB-V452	12s
FS-2 19 mm ¾ circle 	21293328	 Reverse Cutting	Braided	Violet	3-0	30 in (75 cm)	AB-V398	12s
	21293327	 Reverse Cutting	Braided	Violet	4-0	30 in (75 cm)	AB-V397	12s


*Ethilon, Monocryl, PDS II, and Vicryl are trademarks, or registered trademarks of Ethicon Inc.
 **One-Cryl, One-Dox, One-Lon, and One-Sorb are trademarks, or registered trademarks of MWI Veterinary Supply Co.
 ^Caprosyn/Biosyn, Maxon, Monosof, and Polysorb are trademarks, or registered trademarks of Covidien
 ^^Monosorb, Monofil, PGA, and Sorbocryl are trademarks, or registered trademarks of Covetrus
 ^^^Monomend MT, Monomend Max, NY-STA, and Polymend ES are trademarks, or registered trademarks of VPL

Description	Item Code	Needle Type	Material	Suture Color	Size	Length	Mfg Code	Quantity
ECO-MAX™ PDO								
CP-1 36 mm ½ circle 	21293303	▽ Reverse Cutting	Monofilament	Violet	0	30 in (75 cm)	AB-D467	12s
	21293304	▽ Reverse Cutting	Monofilament	Violet	1	30 in (75 cm)	AB-D468	12s
	21293302	▽ Reverse Cutting	Monofilament	Violet	2-0	30 in (75 cm)	AB-D466	12s
CP-2 26 mm ½ circle 	21293307	▽ Reverse Cutting	Monofilament	Violet	0	30 in (75 cm)	AB-D970	12s
	21293306	▽ Reverse Cutting	Monofilament	Violet	2-0	30 in (75 cm)	AB-D969	12s
CT-1 36 mm ½ circle 	21293297	⊙ Taper Point	Monofilament	Violet	0	30 in (75 cm)	AB-D340	12s
	21293296	⊙ Taper Point	Monofilament	Violet	2-0	30 in (75 cm)	AB-D339	12s
CT-2 26 mm ½ circle 	21293295	⊙ Taper Point	Monofilament	Violet	0	30 in (75 cm)	AB-D334	12s
	21293294	⊙ Taper Point	Monofilament	Violet	2-0	30 in (75 cm)	AB-D333	12s
	21293293	⊙ Taper Point	Monofilament	Violet	3-0	30 in (75 cm)	AB-D332	12s
FS-1 24 mm ⅔ circle 	21293300	▽ Reverse Cutting	Monofilament	Violet	2-0	30 in (75 cm)	AB-D451	12s
	21293301	▽ Reverse Cutting	Monofilament	Violet	3-0	30 in (75 cm)	AB-D452	12s
FS-2 19 mm ⅔ circle 	21293299	▽ Reverse Cutting	Monofilament	Violet	3-0	30 in (75 cm)	AB-D398	12s
	21293298	▽ Reverse Cutting	Monofilament	Violet	4-0	30 in (75 cm)	AB-D397	12s
SH 26 mm ½ circle 	21293292	⊙ Taper Point	Monofilament	Violet	2-0	30 in (75 cm)	AB-D317	12s
	21293291	⊙ Taper Point	Monofilament	Violet	3-0	30 in (75 cm)	AB-D316	12s
	21293290	⊙ Taper Point	Monofilament	Violet	4-0	30 in (75 cm)	AB-D315	12s

Description	Item Code	Needle Type	Material	Suture Color	Size	Length	Mfg Code	Quantity
ECO-LON™ NYLON								
FS-26 26 mm $\frac{3}{8}$ circle 	21293179	 Reverse Cutting	Monofilament Nylon	Blue	2-0	30 in (75 cm)	AB-664BL	12s
FS-1 24 mm $\frac{3}{8}$ circle 	21293178	 Reverse Cutting	Monofilament Nylon	Blue	3-0	30 in (75 cm)	AB-663BL	12s
FS-2 19 mm $\frac{3}{8}$ circle 	21293177	 Reverse Cutting	Monofilament Nylon	Blue	4-0	30 in (75 cm)	AB-662BL	12s



Stay Tied in with Aspen

Follow us on Facebook, Instagram and scan our QR code to sign up for e-newsletters so you don't miss a thread of communication from Aspen.



Find us on  

ECO-WEB™

Description	Item Code	Needle Type	Material	Suture Color	Size	Length	Mfg Code	Quantity
ECO-WEB™ PGCL								
CP-1 36 mm ½ circle 	21293316	 Reverse Cutting	Monofilament	Violet	0	36 in (90 cm)	AB-Q967	12s
CT-1 36 mm ½ circle 	21293311	 Taper Point	Monofilament	Violet	2-0	36 in (90 cm)	AB-Q345	12s
	21293310	 Taper Point	Monofilament	Violet	3-0	36 in (90 cm)	AB-Q344	12s
FSL 30 mm ¾ circle 	21293317	 Reverse Cutting	Monofilament	Violet	0	36 in (90 cm)	AB-Q987	12s
FS-1 24 mm ¾ circle 	21293315	 Reverse Cutting	Monofilament	Violet	2-0	36 in (90 cm)	AB-Q943	12s
	21293314	 Reverse Cutting	Monofilament	Violet	3-0	36 in (90 cm)	AB-Q942	12s
FS-2 19 mm ¾ circle 	21293313	 Reverse Cutting	Monofilament	Violet	3-0	36 in (90 cm)	AB-Q923	12s
	21293312	 Reverse Cutting	Monofilament	Violet	4-0	36 in (90 cm)	AB-Q922	12s
SH 26 mm ½ circle 	21293309	 Taper Point	Monofilament	Violet	3-0	30 in (75 cm)	AB-Q316	12s



Aspen Eco-Sutures

EASY AS

1

ORDER

any box of
Aspen
Eco-Suture

2

TRIAL

the sutures with
our risk-free
guarantee!

3

CONSULT

with your
sales specialist
about ordering



Aspen's Eco-Sutures are designed specifically for animal health professionals who value quality and affordability. These sutures are backed by a 100% risk-free guarantee so if you are dissatisfied with the performance in any way, contact your sales representative for a full refund.

Try your first box today!

COMMON USES

Procedure Type	Eco-Suture Choices	Also Possible	SUTURE OPTIONS				
			ANIMAL WEIGHT (lb)				
			15	20	25	30	35
Ligatures	ECO-WEB™, ECO-CRYL™, ECO-MAX™		4-0	4-0 or 3-0		3-0	
Linea Alba	ECO-WEB™ (young) ECO-MAX™ (adult)	ECO-CRYL™		3-0 FS-1, FS-2, CT-2, CT-3, SH, CP-2			← 3-0 or 2-0 →
Intradermal Closure	ECO-WEB™, ECO-CRYL™			4-0 FS-2			← 4-0 or 3-0 →
Subcutaneous Closure	ECO-WEB™, ECO-CRYL™	ECO-MAX™		4-0 SH, FS-2			← 4-0 or 3-0 →
Skin Closure	ECO-LON™		4-0 FS-2		← 4-0 or 3-0 →		3-0 FS-1
Abscesses (Drain)		ECO-LON™				3-0 FS-1	
Anal Saccullectomy	ECO-WEB™	ECO-CRYL™	4-0 FS-2			← 4-0 or 3-0 →	
Biopsies	ECO-LON™	ECO-WEB™, ECO-MAX™	4-0 FS-2			← 4-0 or 3-0 →	
Castration / Cryptorchid	ECO-WEB™, ECO-CRYL™		Cat (none)	3-0 FS-1, FS-2, SH	← 3-0 or 2-0 →		2-0 CP-1, FS-1, FSL
Celiotomy	ECO-WEB™ (young) ECO-MAX™ (adult)	ECO-CRYL™		3-0 FS-1, FS-2, CT-2, SH, CP-2			← 3-0 or 2-0 →
Cesarean Section	ECO-MAX™	ECO-WEB™	4-0 SH	← 4-0 or 3-0 →			3-0 SH, CT-2
Cranial Cruciate Ligament	ECO-MAX™		3-0 SH, CT-2	← 3-0 or 2-0 →			2-0 CP-1, CP-2, CT-1, CT-2, SH
Cystotomy	ECO-WEB™					← 4-0 or 3-0 →	
Dental / Oral Cavity		ECO-WEB™, ECO-CRYL™				← 5-0 or 4-0 →	
Enucleation	ECO-WEB™, ECO-CRYL™					3-0 CT-1, CT-2, CT-3	
Gastric Torsion (gastropexy)	ECO-MAX™			3-0 SH, CT-2			← 3-0 or 2-0 →
Gastrotomy / Enterotomy	ECO-WEB™, ECO-MAX™		4-0 SH			← 4-0 or 3-0 →	
Hernia (umbilical/inguinal)	ECO-MAX™ (recurrent)		3-0 SH, CT-2		← 3-0 or 2-0 →		2-0 CT-1, CT-2, SH
Intestinal Resection / Anastomosis	ECO-WEB™, ECO-MAX™		4-0 SH			← 4-0 or 3-0 →	
Laceration / General Wound	ECO-MAX™ (internal) ECO-CRYL™ (muscular)	ECO-WEB™	4-0 FS-2	← 4-0 or 3-0 →			3-0 FS-2, FS-1
Lumpectomy	ECO-WEB™ (internal) ECO-MAX™ (muscular)		4-0 SH	← 4-0 or 3-0 →			3-0 SH, CT-1, CT-2
Mastectomy	ECO-MAX™		4-0 SH	← 4-0 or 3-0 →			3-0 SH, CT-2
Ovariohysterectomy (OHE)	ECO-WEB™, ECO-CRYL™			3-0 FS-1, FS-2, CT-1, CT-2, CT-3, SH	← 3-0 or 2-0 →		2-0 CT-1, CT-2, FS-1, FSL, CP-1, SH
Splenectomy	ECO-WEB™, ECO-CRYL™		4-0 FS-2			← 4-0 or 3-0 →	

COMMON USES

SUTURE OPTIONS (cont.)									
ANIMAL WEIGHT (lb) (cont.)									
40	45	50	55	60	65	70	75	+	
	3-0		← 3-0 or 2-0 →			2-0	2-0 or 0		0
	2-0 SH, CT-1, CT-2, CP-1, CP-2, FS-1, FSL	← 2-0 or 0 →			0 FSL, CP-1, CP-2, CT, CT-2, CTX	← 0 or 1 →			1 CT-1, CP-1
		3-0 FS-1, FS-2				← 3-0 or 2-0 →			2-0 CP-1, FS-1, FSL
		3-0 SH, CT-2				← 3-0 or 2-0 →			2-0 CT-1, CT-2, SH
	3-0 FS-1	← 3-0 or 2-0 →				2-0 FS			
	3-0 FS-1		← 3-0 or 2-0 →				2-0 FS		
				3-0 FS-1					
				3-0 FS-1					
	2-0 CP-1, FS-1, FSL		← 2-0 or 0 →			0 CP-1		← 0 or 1 →	
	2-0 SH, CT-1, CT-2, FS-1		← 2-0 or 0 →			0 FSL, CP-1, CT-1, CT-2		← 0 or 1 →	
	3-0 SH, CT-2		← 3-0 or 2-0 →			2-0 SH, CT-2, CT-1		← 2-0 or 0 →	
		2-0 CP-1, CP-2, CT-1, CT-2, SH				← 2-0 or 0 →			0 CT-1, CT-2, CP-1, CP-2
	← 4-0 or 3-0 →		3-0 CT-1, SH			← 3-0 or 2-0 →			2-0 CT-1
	← 5-0 or 4-0 →		4-0 FS-2		← 4-0 or 3-0 →		3-0 FS-1, FS-2		← 3-0 or 2-0 →
			3-0 CT-1, CT-2, CT-3						
	← 3-0 or 2-0 →			2-0 CT-1, CT-2, SH					
		3-0 SH, CT-1, CT-2			← 3-0 or 2-0 →				2-0 SH, CT-1, CT-2
	2-0 CT-1, CT-2, SH			← 2-0 or 0 →		0 CT-1, CT-2		← 0 or 1 →	
		3-0 CT-2, SH			← 3-0 or 2-0 →				2-0 SH, CT-1, CT-2
	3-0 FS-2, FS-1		← 3-0 or 2-0 →			2-0 FS-1, FSL			
	3-0 SH, CT-1, CT-2		← 3-0 or 2-0 →			2-0 SH, CT-1, CT-2			
	3-0 SH, CT-2		← 3-0 or 2-0 →			2-0 SH, CT-1, CT-2			
SH	2-0		← 2-0 or 0 →		0 FSL, CP-1, CT-1, CT-2		← 0 or 1 →		1 CT-1, CP-1
		3-0 FS-1, FS-2, SH, CT-2, CT-3				← 3-0 or 2-0 →			2-0 FSL, CP-1, FS-1, CT-1, CT-2, SH



Order through our exclusive distribution partners:

Animal Health 
INTERNATIONAL
a Patterson Company
800.203.5620
animalhealthinternational.com


PATTERSON®
VETERINARY
800.225.7911
pattersonvet.com

Find us on  

05/2022